

**PURITAN MEDICAL PRODUCTS COMPANY LLC****Safety Data Sheet**

Rev. 07 October 27, 2021

**1. Product and company identification****Product name:** Opti-Swab® Liquid Amies Collection & Transport System**Product numbers:** LA-106, LA-116, LA-117**Company identification:**Puritan Medical Products Company LLC  
P.O. Box 149, 31 School Street  
Guilford, Maine 04443-0149 U.S.A.**Contact numbers:**Tel: +1 207-876-3311  
Fax: + 1 207-876-3130**2. Hazards identification**

Skin contact: Irritation

Hazardous ingredients: Irritation

**3. Composition/information on ingredients**

Description:	Liquid Mixture		
Hazardous Components:	Concentration:	R	Classification:
Sodium Chloride	0.3%	36/38	7647-14-5
Calcium Chloride	0.01%%	36	10043-52-4
Potassium Chloride	0.02%%	36/37/38	7447-40-7
Sodium Thioglycolate	0.1%	22/36/37/38	367-51-1

**4. First-aid measures**

Skin contact: Immediately wash with water and soap, then rinse thoroughly.

Eye contact: Flush open eye with plenty of water for 15 minutes. Consult doctor.

Inhalation: Supply fresh air; if breathing is difficult give oxygen. Get medical attention if cough or other symptoms.

Swallowing: Wash mouth with water, if symptoms persist consult doctor.

**5. Fire-fighting measures**

Extinguishing media: Dry chemicals, CO2, foam or Halon extinguisher

Protective Equipment: Personal Protective Equipment (PPE)

**6. Accidental release measures**

Personal precautions: Personal Protective Equipment (PPE)

Environmental precautions: N/A

Methods for cleaning up: Wipe up with damp cloth or mop

**7. Handling and storage**

Handling: No special handling procedures required

Storage: Store away from oxidizing agents. Store at 2°- 25°C.

**8. Exposure controls/personal protection**

Respiratory protection: N/A

Hand protection: Disposable rubber gloves

Eye protection: Safety glasses

Skin and body protection: Protective work clothing (lab coat)

<b>Safety Data Sheet</b>	
<b>9. Physical and chemical properties</b>	
Appearance/Odor	Characteristic
pH:	7.3 ± 0.3 @ 25°C
Density:	Not determined
Boiling point, °C:	Not determined
Melting point, °C	Not determined
Flash point, °C:	Not applicable
Solubility:	N/A
<b>10. Stability and reactivity</b>	
Thermal decomposition/conditions to avoid:	No decomposition and mucous membranes
Incompatible materials:	N/A
Hazardous decomposition products:	No dangerous decomposition products known
<b>11. Toxicology information</b>	
Skin:	Irritant to skin and mucous membranes
Eye:	Irritant effect
Sensitization:	No sensitization effect known
Additional Toxicological Information:	The product shows the following dangers according to internal approved preparations: irritant
<b>12. Ecological Information</b>	
Ecology:	None have been identified.
General Notes:	Not known to be hazardous to water.
<b>13. Disposal considerations</b>	
Recommendation:	Dispose of materials in accordance with federal (40 CFR 261.3), state, and local requirements. This product is not considered a RCRA hazardous waste.
<b>14. Transport information</b>	
ADR/RID	No
IMO/IMDG	No
KAO/IATA:	No
<b>15. Regulatory information</b>	
SARA Section 355:	None of the ingredients is listed
SARA Section 313:	None of the ingredients is listed
TSCA:	7647-14-5 Sodium Chloride 10043-52-4 Calcium Chloride 7447-40-7 Potassium Chloride 367-51-1 Sodium Thioglycolate
California Proposition 65:	None of the ingredients is listed
<b>16. Other information</b>	
Puritan Medical Products Company LLC provides the information in this document in good faith and believes the	

information to be accurate. The chemical, physical and toxicological properties of this product have not been thoroughly investigated. It is the responsibility of the buyer to research and understand safe methods of handling, storing, and disposal of this product. Puritan Medical makes no warranty with respect to such information and assumes no liability for any loss or injury, which may result from the use of this information. It is the buyers responsibility to comply with local, state and federal regulations concerning use and disposal of this product.