25-1503 1PFB SOLID

Puritan® Foam Tipped Applicator

Standard Foam Tip, Solid Polystyrene Handle, Sterile

DESCRIPTION


This specialty applicator is tipped with 100ppi medical grade polyurethane foam. It is bonded securely to a solid polystyrene handle.

FEATURES & BENEFITS

- Medical grade quality, CE marked
- Foam tip has a gentle soft feel creating good collection and release properties
- Can be used dry or with media
- Made in USA
- Product expiration date: 5 years

APPLICATIONS

- Specimen collection for testing and screening
- For nasal and anterior nares specimen collection for diagnostic testing
- Used in Pediatrics,
 Opthalmology, Veterinary Care,
 General purpose use
- Component in rapid diagnostic kits



Item	Tip	Handle	Full	Packaging
No.	length x diameter	length x diameter	Length	pack, weight
non-sterile sterile				
25-1503 1PFB SOLID	16.002 mm x 4.699 mm	73.558 mm x 2.54 mm	76.2 mm	*1,000 BULK
	0.630" x 0.185"	2.896" x 0.100"	3"	4 lbs

Individually wrapped, 1 swab per package - 1,000 Bulk p/case

Products may not appear at actual size. See chart or our website for accurate measurements.

