




Rayon Tipped Applicator- Dry Transport System

Standard tip, polystyrene handle



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kits. This speciality applicator is made with USP grade rayon fiber that is carded on site to maintain consistency. It is bonded securely to a semi-flexible polystyrene or a rigid wood handle by an aqueous based adhesive.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.00" (152.4 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft and smooth tip
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Double swab to accommodate multiple sampling sites
- Made in USA

Applications

- Specimen collection for testing and screening
- Cell collection for DNA testing
- General purpose collection swab



Puritan®

Ref #	Packaging
25-806 2PR BT Sterile	10/50/2 (10 boxes per case of 50 individually wrapped tubes per box, 2 swabs per tube)