



HydraFlock® Flocked Swab - Dry Transport System

Elongated tip, polystyrene handle with breakpoint at 20mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique micro-structure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 5.93" (150.62 mm)

Tip

- 0.219" Dia. (5.56 mm)
- 0.770" Length (19.56 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.89" Length (149.61 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Packaged in a reclosable dry transport tube for ease in transport to lab
- 20mm breakpoint on handle
- Made in USA

Applications

- For collection and release of biological specimens
- Buccal cell collection / DNA testing
- Specimen collection
- Forensic evidence collection
- Rapid diagnostics



Puritan®

Ref #	Packaging
25-3206-H 20MM BT Sterile	10/50/21 (10 boxes per case of 50 individually wrapped tubes per box)