



Environmental Sampling Swab

EnviroMax® Pointed dry foam tip, polypropylene handle



Description

Ready-to-use environmental foam sampling swabs, provided in wet and dry format for accuracy and consistency of sample collection.

Dimensions

6.02" Overall Length (152.908 mm)

Tip

- 0.612" Dia. (15.545 mm)
- 2.220" Length (56.388 mm)

Handle

- 0.125" Dia. (3.175 mm)
- 5.60" Length (142.240 mm)

Qualities

- Large tip for maximized collection and elution
- Quick-turn-open leak-proof cap
- Rigid foam paddle tip
- Large tip for maximized collection and elution
- High-quality Puritan foam-tipped swab
- Shelf-life: 12 months
- Made in USA

Applications

- For environmental surface sampling
- Perfect for spot checks
- Quality control
- Detection of microorganisms

Ref #	Packaging
88010-PF Non-Sterile	25/1 (25 individually wrapped tubes with foam swab)
88010-PFB Non-Sterile	10/25/1 (10 boxes, 25 individually wrapped tubes, 1 swab per package)
25-88010 PF Sterile	25/1 (25 individually wrapped tubes with foam swab)
25-88010 PFB Sterile	10/25/1 (10 boxes, 25 individually wrapped tubes, 1 swab per package)

