



Environmental Sampling Swab

EnviroMax Plus® Round foam tip, Pre-moistened with 50% neutralizing buffer and 50% 0.1% Peptone water, Polypropylene handle



Description

Ready-to-use environmental foam sampling swabs, provided in both wet and dry format for accuracy and consistency of sample collection.

Dimensions

6.02" Overall Length (152.908 mm)

Tip

- 0.612" Dia. (15.545 mm)
- 2.220" Length (56.388 mm)

Handle

- 0.125" Dia. (3.175 mm)
- 5.60" Length (142.240 mm)

Qualities

- Large tip for maximized collection and elution
- Quick-turn-open leak-proof cap
- Swab securely attached to cap
- Rigid foam paddle tip
- High-quality, Puritan foam-tipped swab
- Shelf-life: 12 months
- Made in USA

Applications

- For environmental surface sampling
- Perfect for spot checks
- Quality control
- Detection of microorganisms

Ref #	Packaging
25-88050 PF RND Sterile	25/1 (25 individually wrapped tubes)
25-88050 PFB RND Sterile	10/25/1 (10 boxes, 25 individually wrapped tubes, 1 swab per package)



Puritan®