


SAFETY DATA SHEET

Puritan[®] ESK[®] Environmental Sampling Kit with Neutralizing Buffer 4 mL / 10 mL



 **Puritan Medical Products Co.**
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311
EMAIL sales@puritanmedproducts.com
VISIT puritanmedproducts.com



Puritan® ESK® Environmental Sampling Kit with Neutralizing Buffer 4 mL / 10 mL

1. PRODUCT & COMPANY IDENTIFICATION

Product Name: Puritan® ESK® Environmental Sampling Kit with Neutralizing Buffer, 4 mL / 10 mL

Product Number(s): 25-83004 PD NB, 25-83010 PD NB, 25-83004 PDB NB, 25-83010 PDB NB

Company Identification:

Puritan Medical Products Company
P.O. Box 149, 31 School Street
Guilford, Maine 04443-0149 U.S.A.

Contact numbers:

Tel: +1 207-876-3311
Fax: +1 207-876-3130

2. HAZARDS IDENTIFICATION

This product contains no hazardous constituents, or the concentration of all chemical constituents are below the regulatory threshold limits described by Occupational Safety Health Administration Hazard Communication Standard 29 CFR 1910.1200 and the European Directive 91/155/EEC and 93/112/EC

3. COMPOSITION/INFORMATION ON INGREDIENTS

Description: Liquid mixture

Hazardous Components	Concentration	Classification
N/A	N/A	N/A

4. FIRST-AID MEASURES

Skin contact: Immediately wash with water and soap, then rinse thoroughly

Eye contact: Flush open eye with plenty of water for 15 minutes. Consult doctor.

Inhalation: Supply fresh air; if breathing is difficult give oxygen. Get medical attention if cough or other symptoms appear.

Swallowing: Wash mouth with water, if symptoms persist consult doctor.

5. FIREFIGHTING MEASURES

Suitable extinguishing media: Dry chemicals, CO₂, foam or Halon extinguisher.

Protective Equipment: No special measures required.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures: Not required.

Environmental precautions: Wipe up with damp sponge or mop.

7. HANDLING AND STORAGE

Handling: No special handling procedures required.

Storage: Store away from oxidizing agents.

**8. EXPOSURE CONTROLS/PERSONAL PROTECTION**

Respiratory protection: Not applicable.

Hand protection: Disposable rubber gloves.

Eye protection: Safety glasses.

Skin and body protection: Protective work clothing (Lab coat).

9. PHYSICAL AND CHEMICAL PROPERTIES

Appearance/Odor	Characteristic
pH	7.3± 0.3 @ 25°C
Density	Not determined
Boiling point, °C	Not determined
Melting point, °C	Not determined
Solubility/Miscibility with water	N/A
Flash point, °C	Not applicable

10. STABILITY AND REACTIVITY

Thermal decomposition/conditions to avoid: No decomposition if used according to specifications.

Incompatible materials: N/A

Hazardous decomposition products: N/A

11. TOXICOLOGY INFORMATION

Skin: No irritant effect.

Eye: No irritant effect.

Sensitization: No sensitization effects known.

Additional Toxicological Information

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

12. ECOLOGICAL INFORMATION

Ecotoxicological Effects: None have been identified.

General Notes: Not known to be hazardous to water.

13. DISPOSAL CONSIDERATIONS

Recommendation: Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements. This product is not considered a RCRA hazardous waste.

14. TRANSPORT INFORMATION

ADR/RID: No.

IMO/IMDG: No.

KAO/IATA class: No.

**15. REGULATORY INFORMATION**

SARA Section 355: None of the ingredients are listed.

SARA Section 313: None of the ingredients are listed.

TSCA: N/A

California Proposition 65: None of the ingredients are listed.

16. OTHER INFORMATION

Puritan Medical Products Company LLC provides the information in this document in good faith and believes the information to be accurate. The chemical, physical and toxicological properties of this product have not been thoroughly investigated. It is the responsibility of the buyer to research and understand safe methods of handling, storing, and disposal of this product. Puritan Medical makes no warranty with respect to such information and assumes no liability for any loss or injury, which may result from the use of this information. It is the buyers responsibility to comply with local, state and federal regulations concerning use and disposal of this product.



Puritan®
Quality since 1919



Puritan Medical Products Co.
31 School Street, P.O. Box 149
Guilford, Maine 04443 USA

CALL 800.321.2313 • 207.876.3311

EMAIL sales@puritanmedproducts.com

VISIT puritanmedproducts.com