



ESK - Letheen Broth

ESK® Sampling Kit - Pre-filled with 4 ml or 10 ml Letheen Broth with Tween 80



Description

Letheen Broth is recommended for testing materials for preservatives where neutralization of disinfectants such as quaternary ammonium is important for antimicrobial activity.

Dimensions

4.06" Overall Length of Kit (103.124 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Handle

- 3.25" Length (82.55 mm)

Qualities

- Easy edge-peel packaging
- Quick-turn leak-resistant cap
- High-quality Puritan polyester tipped swab
- Swab is securely attached to cap
- Shelf-life: 18 months
- Made in USA

Applications

- For environmental surface sampling
- Quality Control
- Detection of Microorganisms

Ref #	Packaging
25-83004 PD LB Sterile	50/1 (4 ml - 50 individually wrapped kits with polyester swab)
25-83004 PDB LB Sterile	6/50/1 (4 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)
25-83010 PD LB Sterile	50/1 (10 ml - 50 individually wrapped kits with polyester swab)
25-83010 PDB LB Sterile	6/50/1 (10 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)



Puritan®