



ESK - Neutralizing Buffer

ESK® Sampling Kit - Pre-filled with 4 ml or 10 ml Neutralizing Buffer



Description

An ideal choice for on-the-spot quality control checks. The Neutralizing Buffer is recommended for verification in your sanitation programs with equipment, according to the standard methods.

Dimensions

4.06" Overall Length of Kit (103.124 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.450 mm)

Handle

- 3.25" Length (82.55 mm)

Qualities

- Easy edge-peel packaging
- Quick-turn leak-resistant cap
- High-quality Puritan polyester-tipped swab
- Swab is securely attached to cap
- Shelf-life: 18 months
- Made in USA

Applications

- For environmental surface sampling
- Quality control
- Detection of microorganisms

Ref #	Packaging
25-83004 PD NB Sterile	50/1 (4 ml - 50 individually wrapped kits with polyester swab)
25-83004 PDB NB Sterile	6/50/1 (4 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)
25-83010 PD NB Sterile	50/1 (10 ml - 50 individually wrapped kits with polyester swab)
25-83010 PDB NB Sterile	6/50/1 (10 ml - 6 boxes, 50 individually wrapped kits with polyester swab, 1 swab per package)