

PURITAN MEDICAL PRODUCTS COMPANY LLC

Safety Data Sheet

Rev. 05 November 17, 2020

1. Product and company identification	
Product Name: Neutralizing Buffer 4mL	
Product Number: 25-83004 PD NB	
Company Identification: Puritan Medical Products Company LLC P.O. Box 149, 31 School Street Guilford, Maine 04443-0149 U.S.A.	Contact numbers: Tel: +1 207-876-3311 Fax: + 1 207-876-3130
2. Hazards identification	
This product contains no hazardous constituents, or the concentration of all chemical constituents are below the regulatory threshold limits described by Occupational Safety Health Administration Hazard Communication Standard 29 CFR 1910.1200 and the European Directive 91/155/EEC and 93/112/EC.	
3. Composition/information on ingredients	
Description:	Liquid Mixture
Hazardous Components:	N/A
Concentration:	N/A
Classification:	N/A
4. First-aid measures	
Skin contact:	Immediately wash with water and soap, then rinse thoroughly.
Eye contact:	Flush open eye with plenty of water for 15 minutes. Consult doctor.
Inhalation:	Supply fresh air; if breathing is difficult give oxygen. Get medical attention if cough or other symptoms appear.
Swallowing:	Wash mouth with water, if symptoms persist consult doctor.
5. Fire-fighting measures	
Suitable Extinguishing media:	Dry chemicals, CO ₂ , foam or Halon extinguisher.
Protective Equipment	No special measures required.
6. Accidental release measures	
Personal precautions, protective equipment and emergency procedures:	Not required
Environmental precautions:	Wipe up with damp sponge or mop.
7. Handling and storage	
Handling:	No special handling procedures required
Storage:	Store away from oxidizing agents.
8. Exposure controls/personal protection	
Respiratory protection:	Not Applicable
Hand protection:	Disposable rubber gloves
Eye protection:	Safety glasses
Skin and body protection:	Protective work clothing (Lab coat)

9. Physical and chemical properties	
Appearance/Odor:	Characteristic
pH:	7.3± 0.3 @ 25°C
Density:	Not determined
Boiling point, °C:	Not determined
Melting point, °C	Not determined
Solubility/Miscibility with water:	N/A
Flash point, °C:	Not applicable
10. Stability and reactivity	
Thermal decomposition/conditions to avoid:	No decomposition if used according to specifications
Incompatible materials:	N/A
Hazardous decomposition products:	N/A
11. Toxicology information	
Skin:	No irritant effect
Eye:	No irritant effect
Sensitization	No sensitization effects known
Additional Toxicological Information:	When used and handled according to specifications, the product does not have any harmful effects according to our
12. Ecological Information	
Ecotoxicological Effects:	None have been identified
General Notes:	Not known to be hazardous to water
13. Disposal considerations	
Recommendation:	Dispose of material in accordance with federal (40 CFR 261.3), state and local requirements. This product is not considered a RCRA hazardous waste.
14. Transport information	
ADR/RID:	No
IMO/IMDG:	No
KAO/IATA class	No
15. Regulatory information	
SARA Section 355:	None of the ingredients are listed.
SARA Section 313:	None of the ingredients are listed.
TSCA:	N/A
California Proposition 65:	None of the ingredients are listed.
16. Other information	
Puritan Medical Products Company LLC provides the information in this document in good faith and believes the information to be accurate. The chemical, physical and toxicological properties of this product have not been thoroughly investigated. It is the responsibility of the buyer to research and understand safe methods of handling, storing, and disposal of this product. Puritan Medical makes no warranty with respect to such information and assumes no liability for any loss or injury, which may result from the use of this information. It is the buyers responsibility to comply with local, state and federal regulations concerning use and disposal of this product.	