





















Packaging Symbol Glossary

Symbol	Meaning
	Catalogue number
	Batch Code
	Use by date
	Caution, consult accompanying documents
	Do not re-use, single use
	Date of manufacture
	Manufacturer
	CE mark. Product conforms with the essential requirements in the European Medical Devices Directive 93/42/EEC.
	CE mark including SGS (notified body) identification number. Product conforms with the essential requirements in the European Medical Devices Directive 93/42/EEC.
	Product does not contain latex

Symbol	Meaning
	Sterilized using Ethylene Oxide
	Sterilized using irradiation
	Do not use if package is damaged
	Non-sterile
	Authorized Representative in the European Community
	Temperature limitation
	Consult Instructions for Use
	In vitro diagnostic medical device
	Caution, warning
	Quantity of devices