

Puritan	RELATED PRODUCTS					PRODUCT DESCRIPTION			
	Copan	BD	Hardy Diagnostics	Medical Wire	Remel: ThermoFisher Sci	Description	Vial Diameter x Length (mm)	Fill Volume	Cap Color
BHI-300	474CE.A	221813 (5mL)	R15 (2mL; 13 x 100 mm tube)	N/A	R060264/R060266 (4mL)	BHI Broth: Enrichment medium that facilitates growth of fastidious aerobic bacteria. For In Vitro Diagnostic Use.	12 x 80	2mL	purple
GN-200	085CU.A (4mL)	221730 (8mL)	K01 (6mL; 15 x 103 mm)	N/A	R060990/R060992 (5mL)	GN Broth: Enrichment broth medium used for the isolation of <i>Salmonella</i> spp and <i>Shigella</i> spp. For In Vitro Diagnostic Use.	12 x 80	2mL	black
TSB-200	4E023N.A	221351	K126 (5mL)	N/A	R065032/R065030 (5mL)	TSB Salt Broth: Enrichment broth medium used for the isolation of specimens containing MRSA. For In Vitro Diagnostic Use.	12 x 80	2mL	blue
CB-200	n/a	n/a	n/a	n/a	n/a	Cary-Blair: Transport medium for the preservation of fecal and rectal specimens containing enteric pathogenic bacteria. For In Vitro Diagnostic Use.	12 x 80	2mL	orange
CB-500	n/a	n/a	n/a	n/a	R060450/ R060452 (9mL)		16 x 100	5mL	
CT-200	094CE.A	n/a	n/a	n/a	n/a	CT Broth: Selective enrichment broth medium for use in the isolation of <i>Candida</i> spp and <i>Trichomonas</i> spp. For in vitro diagnostic use.	12 x 80	2mL	blue
CT-500	092CN (4ml)	298323	n/a	n/a	n/a		16 x 100	5mL	
LB-200	476CE.A	n/a	n/a	MWGBST	R064812 (3mL)	Lim Broth: Selective enrichment medium for use in the qualitative procedures for the isolation of Group B Streptococcus (GBS) from clinical specimens. For In Vitro Diagnostic Use.	12 x 80	2mL	light blue
LB-500	n/a	292209/ 296266	L57	n/a	R064814		16 x 100	5mL	
MRSA-200	477CE.A (TSB Salt Broth)	n/a	n/a	n/a	n/a	MRSA: Enrichment broth medium used for the isolation of specimens containing MRSA. For In Vitro Diagnostic Use.	12 x 80	2mL	pink
MRSA-500	n/a	n/a	n/a	n/a	n/a		16 x 100	5mL	

RL-200	n/a	n/a	n/a	n/a	n/a	Regan-Lowe: Semi-solid transport medium for the isolation of <i>Bordetella</i> species. For In Vitro Diagnostic Use.	12 x 80	2mL	white
RL-500	n/a	297855 (10mL)	Q32 (4mL; 13 x 100mm tube)	n/a	R064141		16 x 100	5mL	
SB-200	475CE.A	n/a	n/a	MWSEL	n/a	Selenite Broth: Enrichment broth medium used for the isolation of <i>Salmonella spp.</i> and <i>Shigella spp.</i> For In Vitro Diagnostic Use.	12 x 80	2mL	white
SB-500	n/a	221020 (8mL)	n/a	n/a	R064500 / R064502 (9mL)		16 x 100	5mL	
SP-100	097CE.A	n/a	n/a	n/a	n/a	Sputum Solution: Liquefying solution of clinical sputum specimens for the isolation of bacteria and fungi that cause respiratory tract infections. For In Vitro Diagnostic Use.	12 x 80	1mL	natural
SP-100P	098CE.A	n/a	n/a	n/a	n/a		12 x 80	1mL	
VL 1201 SAFE	608C	n/a	n/a	n/a	n/a	PurSafe® DNA/RNA Preservative: Transport medium for the preservation of DNA and RNA of clinical specimens. Not for diagnostic use. Research use only.	12 x 80	1mL	light blue
VL 1201 SAFE-H	608CS01R	n/a	n/a	n/a	n/a		12 x 80	1mL	
VL 1202 SAFE	606C	n/a	n/a	n/a	n/a		12 x 80	2mL	
VL 1202 SAFE-H	606CS01R	n/a	n/a	n/a	n/a		12 x 80	2mL	