



Puritan®

Quality since 1919





Puritan[®] ENVIRONMENTAL

This line is a direct result of collaborating with clients in food production, pharmaceuticals, and cosmetics. Puritan[®] environmental kits meet the needs of industries that require ongoing sampling to meet rigorous quality control standards.

Puritan's EnviroMax[®], EnviroMax Plus[®] and ESK[®] Sampling Kits provide all-in-one ease of use. They offer oversized, paddle-type swabs for flat surfaces as well as smaller, shaped tips to reach into corners and crevices that might otherwise be overlooked.

All the swabs in our environmental line have white handles and tips and are packaged in a translucent plastic vial.

EnviroMax® & EnviroMax Plus®

The ultimate in environmental surface sampling from the swab experts at Puritan.


The oversized foam swabs are strong but flexible – perfect for large surfaces and equipment parts.



IMPROVED DESIGN!

LEAK-PROOF CAP

No loss of sample – both tube and cap are leak resistant to block gasification and prevent liquid spills.



PACKAGING OPTIONS

Available both individually wrapped and unwrapped.

STRONG HANDLE

Allows controlled pressure to the paddle-like foam swab.

MAXIMIZED COLLECTION

Absorbent foam tips available pre-moistened or dry – your choice depending on your application.

POINTED OR ROUNDED

There's a swab for every application.





FEATURED PRODUCT:
EnviroMax® Sampling Kit

puritanmedproducts.com

Products may not appear at actual size. See chart dimensions or our website for accurate measurements. 3/2019





ENVIROMAX® SAMPLING SWABS (DRY)

Produced as a complete surface-testing kit. Pointed tip for crevices and corners.
Rounded tip for counters, floors, and larger equipment surfaces.

	Item No.	Wrapping Style	Product Description	Tip length x diameter	Swab Length	Packaging pack, weight
non-sterile sterile						
 	88010-PF	wrapped	Pointed foam tip, polypropylene handle	56.388 x 15.545 mm 2.22 x 0.612"	147.32 mm 5.8"	25/1, 2 lbs
	25-88010 PF	wrapped				
	25-88010 PF UW	unwrapped				
	88010-PFB	wrapped				10/25/1, 20 lbs
	25-88010 PFB	wrapped				
	25-88010 PFB UW	unwrapped				
 	88010-PF RND	wrapped	Rounded foam tip, polypropylene handle	56.388 x 15.545 mm 2.22 x 0.612"	145.034 mm 5.71"	25/1, 2 lbs
	25-88010 PF RND					
	88010-PFB RND					10/25/1, 20 lbs
	25-88010 PFB RND					

ENVIROMAX PLUS® SAMPLING SWABS (PRE-MOISTENED)

Produced as a complete surface-testing kit. The pre-moistened swab is integrated with the handle and cap for ease and accuracy of sampling. Pre-moistened with 50% Neutralizing Buffer and 50% 0.1% Peptone Water.

	Item No.	Wrapping Style	Product Description	Tip length x diameter	Swab Length	Packaging pack, weight
non-sterile sterile						
 	25-88050 PF	wrapped	Pointed foam tip, polypropylene handle	56.388 x 15.545 mm 2.22 x 0.612"	147.32 mm 5.8"	25/1, 2 lbs
	25-88050 PF UW	unwrapped				
	25-88050 PFB	wrapped				10/25/1, 20 lbs
	25-88050 PFB UW	unwrapped				
 	25-88050 PF RND	wrapped	Rounded foam tip, polypropylene handle	56.388 x 15.545 mm 2.22 x 0.612"	145.034 mm 5.71"	25/1, 2 lbs
	25-88050 PFB RND	wrapped				10/25/1, 20 lbs


Did you know?


The CDC's **Vessel Sanitation Program (VSP)** includes an inspection component for cruise ships. Periodic, unannounced operational sanitation inspections incorporate testing for the presence of norovirus to control the introduction, transmission, and spread of this gastrointestinal illness.

ENVIROMAX® AND ENVIROMAX PLUS® NV SAMPLING SWABS

Used for the detection of human norovirus on environmental surfaces, utilizing a newly developed protocol by the Centers for Disease Control to assist in determining the level of environmental contamination during outbreaks.

FEATURED PRODUCT: *EnviroMax® NV Sampling Swab (#25-88060 PF UW DRY)*



	Item No.	Wrapping Style	Product Description	Tip length x diameter	Swab Length	Packaging pack, weight
non-sterile sterile						
	25-88060 PF UW	unwrapped	Pre-moistened, large round foam tip, polypropylene handle	34.925 x 21.59 mm 1.375 x 0.850"	147.32 mm 5.8"	25/1, 2 lbs
	25-88060 PF UW DRY	unwrapped	Dry, large round foam tip, polypropylene handle			


Did you know?

The CDC estimates that each year **roughly 1 in 6 Americans (or 48 million people) gets sick**, 128,000 are hospitalized, and 3,000 die of foodborne diseases.




WHITE HANDLE & TIP, TRANSLUCENT VIAL

All the swabs in our environmental line have white handles and tips and are packaged in a translucent plastic vial. Reactions – even small, early reactions – are visible, saving you time and effort.



CATCH & RELEASE

Foam tips capture an ample volume of material and release readily.




PADDLE STRUCTURE

Oversized tips with a rigid paddle-like structure within the foam gives the tip stability and makes it easy to reach into crevices where foam-only tips may bunch up, crush, or fold over.

ESK® Sampling Kits

You can count on accurate, repeatable laboratory results.

Simple, easy-to-use sampling kits are recommended for environmental monitoring, hygiene control, sanitation programs, and bioburden determination. Each sterile ESK® comes in a frustration-free peel pouch.



FEATURED PRODUCT:

ESK® Sampling Kit

puritanmedproducts.com

Products may not appear at actual size. See chart dimensions or our website for accurate measurements. 3/2019


ESK® SAMPLING KITS

The swab is attached to the cap, giving you extra, rigid handle length for firm control. When finished, simply put the swab back in the tube and close the cap tightly. Choose from 4 standard options (in either 4 ml or 10 ml configurations) or request custom media or fill volumes.

	Item No.	Product Description	Tip length x diameter	Vial diameter x length	Fill Volume, Medium	Packaging pack, weight
	non-sterile sterile					
	25-83004 PD BPW	Polyester tip, polystyrene handle	17.45 x 5.156 mm 0.687 x 0.203"	16 x 100 mm	4 ml, Buffered Peptone Water	50/1, 3 lbs
	25-83004 PDB BPW					6/50/1, 18 lbs
	25-83010 PD BPW				10 ml, Buffered Peptone Water	50/1, 3 lbs
	25-83010 PDB BPW					6/50/1, 19 lbs
	25-83004 PD BS	Polyester tip, polystyrene handle	17.45 x 5.156 mm 0.687 x 0.203"	16 x 100 mm	4 ml, Butterfield's Solution	50/1, 3 lbs
	25-83004 PDB BS					6/50/1, 18 lbs
	25-83010 PD BS				10 ml, Butterfield's Solution	50/1, 3 lbs
	25-83010 PDB BS					6/50/1, 19 lbs
	25-83004 PD LB	Polyester tip, polystyrene handle	17.45 x 5.156 mm 0.687 x 0.203"	16 x 100 mm	4 ml, Lethen Broth	50/1, 3 lbs
	25-83004 PDB LB					6/50/1, 18 lbs
	25-83010 PD LB				10 ml, Lethen Broth	50/1, 3 lbs
	25-83010 PDB LB					6/50/1, 19 lbs
	25-83004 PD NB	Polyester tip, polystyrene handle	17.45 x 5.156 mm 0.687 x 0.203"	16 x 100 mm	4 ml, Neutralizing Buffer	50/1, 3 lbs
	25-83004 PDB NB					6/50/1, 18 lbs
	25-83010 PD NB				10 ml, Neutralizing Buffer	50/1, 3 lbs
	25-83010 PDB NB					6/50/1, 19 lbs

Did you know?

Our ESK® product line delivers an **easy-to-use kit** that allows for sampling from contact points and tight spaces in every area of a production facility, with minimal interruption of operations.



TIGHT SPACE FIT

The swab is securely attached to the cap for confident control and ease of use. Ideal for collecting samples from corners, cracks and crevices.

Puritan® – a tradition of excellence

Order with confidence, knowing you will always get consistent quality, ready availability, and outstanding service. You can count on Puritan.

From our inception in 1919, our wood products operation has grown to become a state-of-the-art medical device manufacturer serving markets worldwide with Puritan's US-made products.

Our success is evident in the success of our customers. Our guarantee is to produce the highest-quality single-use products possible. Our mission is to provide our customers with unparalleled service. Our promise is to deliver – the right product, on time.

- **CALL**

800-321-2313 (US & Canada)
207-876-3311

- **VISIT**

31 School Street
PO Box 149
Guilford, Maine 04443

- **CONNECT**

puritanmedproducts.com

- **EMAIL**

sales@puritanmedproducts.com

Puritan Medical Products Co.
31 School Street
PO Box 149
Guilford, Maine 04443-0149 USA

800-321-2313 (US & Canada)
207-876-3311
puritanmedproducts.com
sales@puritanmedproducts.com