



# Clean Earth Technologies, LLC

101 North Chestnut Street, Suite 101

Winston-Salem, NC 27101

336-397-1900 (T); 336-397-1911 (F)

Form# LS1000

## Certificate of Analysis for Forensic Swabs

Test Material: Foam-tipped Swab

Batch and lot details: PSR# 207-09, Lot 2198

Storage: Ambient

Test Dates: October 9, 2009 to October 14, 2009

Test Description	Method	Results*
DNase contamination of swabs	Contaminant enzyme was extracted from swabs, incubated with DNA and analyzed by gel electrophoresis. The relative intensity of negative control and sample bands must be equal.	Pass- No contamination
Human DNA contamination of swabs	Contaminant human DNA was extracted from swabs, specific DNA segment was PCR amplified, and analyzed by gel electrophoresis.	Pass- No contamination

\*under the test conditions

Tested by: Julie Turner

Date: 10-16-09

Approval: [Signature]

Date: 10/16/2009

**CET, LLC Proprietary**

Issued Oct 15, 2009